

CHCP-I-3D Cynomolgus Hepatocytes Cryopreserved Plateable for Induction assays and for 3D culture

Cell Specification – Certificate of Analysis (CoA)

Lot CH240516 Batch Release: July 29, 2024

Donor data

Species: Macaca fascicularis Gender: female

Age: 7 years 3 months

Cryopreservation and Thawing

Cryopreservation:

Date: May 16, 2024

Amount per vial: 10.0 x 10⁶ cells

Thawing: n=1


Post-thaw viability: 95.9 %

Post-thaw yield per vial: 6.5 x 10⁶ cells


Recovery: 65 %

2D culture

Phase contrast on day 1 after thawing (24well plate)



Phase contrast on day 3 after thawing (24well plate)



Recommended seeding density on collagen-coated plates:

24well plate - 400,000 cells/well // 96well plate - 80,000 cells/well.

Culture in Human Hepatocyte Maintenance Medium (HHMM).



CYP P450 activity in 2D culture after thawing:

Ethoxyresorufin-O-deethylation:

Induction with 25 μM β-Naphthoflavone

 $pmol/(mg \times min)$ 24well: 42.9 ± 3.6

96well: 89.9 ± 25.5

X-fold induction

3.4


4.7

3D culture


Cells seeded in 96well ULA round bottom plates (FaCellitate), 2,500 cells/well

day 8

day 10



scale bar 100 µm



scale bar 100 µm

Suspension culture

Viability test on orbital shaker (Eppendorf Thermomixer C, 1000 rpm at 37 °C with 0.5 x 10⁶ cells in 0.5 ml HPM-Cryo):

| Time (h) | 0 | 1 | 2 | 3 | 4 | 5 |
|---------------|------|------|------|------|------|------|
| Viability (%) | 95.9 | 89.6 | 89.3 | 91.0 | 89.2 | 87.6 |

Note: Yield, viability, recovery and activity assays were performed at PRIMACYT using PRIMACYT's manual for thawing, plating and culture of primary cryopreserved hepatocytes.

Store at -150 °C or in the vapour phase of LN₂

This product should be considered as potential biohazard. Only intended for *in vitro* use.

| sued by: A. Sattler | Verified by: T. Krimmling |
|---------------------|---------------------------|
|---------------------|---------------------------|