

Rainbow Trout Liver Subcellular fractions Product specification – Certificate of Analysis (CoA)					
Lot RTL230209 (P	ool of 10)	Batch Release: May 15, 2025			
Product information					
Product number	Product descrir	Product description			

Product number	Product description	Amount	Protein content
RTL-S9-2P10	Rainbow trout liver S9 fraction, female, pool of 10	0.5 mL	20 mg/mL

Donor data				
Species: Rainbow trout (Oncorhynchus mykiss)	Gender: female Age: sexual immature Pool: N = 10			

Equal amounts of Liver tissues were pooled to generate subcellular fractions.


Animal characteristics										
Donor	1	2	3	4	5	6	7	8	9	10
Fish weight (g)	550	537	531	642	460	651	546	439	424	588
Fish length (cm)	34	34	36	37	33	38	35	34	33	36
Gonaden weight (g)	1.4	1.1	2.0	0.8	0.5	1.8	0.7	1.4	0.4	0.8
GSI (gonaden weight/fish weight)	0.25	0.20	0.38	0.12	0.11	0.28	0.13	0.32	0.09	0.14
liver weight (g)	72.2									

GSI = Gonadosomatic index



Enzyme assay results					
Enzyme (Human isoforms)	Assay	Enzyme activities (nM/min)			
	Assay	S9 Fraction			
CYP1A2	Phenacetin-O-deethylase	0.74 ± 0.17			
CYP2A6	Coumarin-7'-hydroxylase not detectable				
CYP2B6	Bupropion-hydroxylase	0.08 ± 0.02			
CYP2C9	Diclofenac 4'-hydroxylase	3.51 ± 0.31			
CYP2C19	Mephenytoin 4'-hydroxylase	not detectable			
CYP2E1	Chlorzoxazone 6'-hydroxylase	0.69 ± 0.13			
CYP3A4	Midazolam 1'-hydroxylase	0.36 ± 0.11			
UDP-GT	UDP-Glucuronosyltransferase	313.38 ± 75.89			
SULT	Sulfotransferase	10.54 ± 1.81			

Note: Activity assays were performed at PRIMACYT GmbH. The assays were conducted at 1 mg/mL protein in 0.1 M Phosphate buffer at 14 °C for 15 min (phase I) and 30 min (phase II). Values are expressed as mean ± SD of 2 separate experiments.



Note: SDS-PAGE was performed at PRIMACYT GmbH using a 7 % stacking gel and a 10 % separating gel with 10, 20 and 40 µg protein as indicated. Pageruler Prestained Protein Ladder by Thermo Scientific was used as marker. Protein sizes (kDa) denoted on the left.

Store at -80 °C.

This product should be considered as potential biohazard. Only intended for *in vitro* use.

Issued by: M. Reu

Verified by: T. Krimmling

© PRIMACYT Cell Culture Technology GmbH, Hagenower Str. 73, D-19061 Schwerin, FRG, www.primacyt.com E-mail: info@primacyt.com, Tel.: +49 (0) 385-3993-600, Fax: +49 (0) 385-3993-602